



PINBALL GAME FOAM CUSHION PACKING


▼ Left Side Foam Cushions ▼



▼ Right Side Foam Cushions ▼



- 3: BACKBOX FRONT TOP RIGHT
- 4: BACKBOX FRONT TOP LEFT
- 1: BACKBOX FRONT BOTTOM LEFT
- 2: BACKBOX FRONT BOTTOM RIGHT
- 10: CABINET FRONT BOTTOM RIGHT
- 9: CABINET FRONT BOTTOM LEFT




◀ Front Right Views ▼

Position # 1 as shown under Right Side Bottom of Backbox before lowering all the way



Front Left View ▶

Position # 2 as shown under Left Side Bottom of Backbox before lowering all the way




Save all packing materials for future use.

When transporting your game, re-pack with all foam cushions. Each foam cushion is numbered and can only fit **properly** one way; *odd numbers on the right and even numbers on the left* (view the photos for correct placement).

Ensure the four (4) short cardboard "V-Guards" and the two (2) long cardboard "V-Guards" are positioned under the strapping to protect the edges of your game.

▲ Rear Left Views ▲

- 8: CABINET REAR TOP LEFT
- 7: CABINET REAR TOP RIGHT
- 6: CABINET REAR BOTTOM LEFT
- 5: CABINET REAR BOTTOM RIGHT



In the shipping box, the game legs are stored in two (2) square tubes along each bottom edge of the pinball game cabinet.

1-800-542-5377



Stern® Pinball, Inc. © 2004
SPI Part № 755-5310-01